

HUBBARD COMMUNICATIONS OFFICE
Saint Hill Manor, East Grinstead, Sussex

HCO BULLETIN OF 11 FEBRUARY 1966

Remimeo
All Students
All Sen Staff
Franchise

FREE NEEDLES, HOW TO
GET THEM ON A PC

Free needles can be obscured only by overruns and auditor goofs in the rehab session and ARC Breaks in past auditing.

When a TA goes up or is up it means an overrun in life or on a process or grade of release.

The only place you can't get an overrun is at Grade VII. All grades below that are subject to overrun.

Life subjects are subject to overrun before Scientology. The mechanism is this: one conceived a purpose. He or she succeeded in it, then kept on and overran it. In auditing one hits the purpose and the overrun of it and gets a free needle on it. That doesn't mean the person was a release then. It means that the spotting of the purpose and the overrun by auditing produces a free needle today.

It may be necessary to find whole track overruns on some pcs in rehabilitation of grades. If a lot of levels have been run past free needle it may be necessary to take apart the mess like a bundle of yarn to get the first free needle. In such a case one rehabs any grade the pc has been run on that the pc can remember. One handles this briefly until the pc is happy but not necessarily to free needle. One then finds another overrun, does the same. One goes on and on looking for moments the pc felt good about processing at one or another time. If you keep this up, suddenly you will see a free needle on the pc! Establish what grade it is free on, then quickly get the needle free on the remaining overrun grades (but not grades pc was never run on). It may be necessary to take into account a whole track overrun of a purpose or even the purpose to get release, clear or O.T.

It is all very quick, deft auditing, very much on procedure using standard rehab tech - but no repetitive grind.

You won't see a freeing up of a needle unless you set your sensitivity on a Mark V to a stiff needle for the pc. You can increase sensitivity or decrease it as the pc progresses, but by setting the sensitivity so the needle is pretty still and stiff you will see easily a freeing up of the needle and then a free needle. Using sensitivity T28 will obscure every free needle as the needle is too loose already for the auditor to see any change.

Pcs are most apt to go free needle after a big cog. So don't be so engrossed in looking at the pc during cognitions. Keep an eye on that needle. And if it goes free, don't ask anything else. Just gently give the pc a "That's it" and without a chop of comm, ease the pc off to "Declare?" in Qual. (Or if a field auditor, start the next grade.)

Gently, gently, smooth TRs get you free needles.

A dirty needle is always caused by auditor chops, flubs etc. You can always trace a dirty needle right back to a TR error by the auditor. If a needle goes dirty in a rehab session, get the List I out right now and quickly find why. It's always an auditor goof on the TRs or Tech procedure.

Rehabs are not a substitute for processes. If a grade hasn't been run, you can't rehab it of course.

In rehab, never use a new process to cure an overrun. Rehab the process that was overrun, not new ruds.

And see HCO Pol Ltr 10 Feb 1966 on this subject.

You can get free needles on pcs. It just requires standard TRs, standard tech, standard rehab and wanting to get one and letting a pc have one.